

Apostle MiniGenomics® Whole Blood Genomic DNA Isolation Kit, Automation Protocol

Isolation of genomic DNA from Whole Blood (350 µL) with MagTouch 2000, 96-well

Revision A.1

Product Description

The Apostle MiniGenomics® Whole Blood Genomic DNA Isolation Kit is designed for rapid isolation of high-quality, ready-to-use genomic DNA from fresh or frozen whole blood containing Citrate, EDTA. The kit uses proprietary Apostle MiniGenomics® technology and offers highly efficient, reproducible recovery of high-quality DNA with high yield. This kit is designed for input volumes ranging from 50 µL to 350 µL of whole blood. The isolated DNA samples are suitable for a broad range of subsequent applications, including sequencing, PCR, etc. The protocol is designed for 96-well plate automated on Apostle MagTouch 2000 Automation Platform.

Required materials

Apostle MiniGenomics® Whole Blood Genomic DNA Isolation Kit (Cat# A230814-50 or A230814-200)

Apostle MagTouch 2000 (Cat# A201126-96)

Apostle 96-well deep well plates and tip comb

DNase/RNase-free tubes (1.5 mL and 15 mL)

Adjustable single and multi-channel micropipettes and tips (1000 µL, 200 µL, and 20 µL)

Automated

1. Set up the 96-well plates according to the table below, outside the instrument
2. Select the [WBLy350] program. Open the cabin door, and place the plates into the MagTouch 2000 as indicated on the instrument display. Press the position button to turn the rotary table and place all the plates in turn.
3. After placing all the plates, press the [Run] icon on the interface to start the program.

Plate Position	Plate Type	Name	Content	Reagent volume for each well
1	96 well plate with comb	-Load-	Comb	-
2	96 well plate	Lysis	Proteinase K	35 µL
			Whole blood sample	350 µL
			Lysis Enhancer Solution	20 µL
			Lysis/Binding Solution	170 µL

Vortex or shaker

Tabletop centrifuge

Heater (for sample lysis)

Ethanol, 200 proof

80% ethanol

Nuclease-free water

Procedure for isolation of gDNA

A. Sample lysis

Manual

1. Add 35 µL Proteinase K, vortex briefly and centrifuge briefly.
2. Add 350 µL of whole blood sample to the tube, vortex for 1 to 3 mins, and centrifuge briefly.
3. Add 20 µL Lysis Enhancer solution to the tube.
4. Add 170 µL of Lysis Binding solution, vortex and centrifuge briefly.
5. Incubate the tube at 60°C for 10 mins.

B. DNA extraction (automated)

1. **(Optional)** Mix Magnetic Nanoparticles in ethanol, 200 proof with a 1:76 ratio before use. The amount of binding mixture required is 385 μL per sample.
2. Set up the 96-well plates according to the table below, outside the instrument

Plate Position	Plate Type	Name	Content	Reagent volume for each well
1	96 well plate with comb	-Load-	Comb (same as in previous step)	-
3	96 well plate	Bind	Lysed sample from step A	575 μL
			Ethanol, 200 proof	380 μL
			Magnetic Nanoparticles	5 μL
			Binding Enhancer	4 μL
4	96 well plate	Wash1a	Wash Solution	1000 μL
5	96 well plate	Wash1b	Wash Solution	1000 μL
6	96 well plate	Wash2a	80% Ethanol, molecular biology grade	1000 μL
7	96 well plate	Wash2b	80% Ethanol, molecular biology grade	1000 μL
8	96 well plate	Elu	Elution Solution	100 μL

3. Select the [WBEx350] program. Open the cabin door, and place the plates into the MagTouch 2000 as indicated on the instrument display. Press the position button to turn the rotary table and place all the plates in turn.
4. After placing all the plates, press the [Run] icon on the interface to start the program.
5. When the program is complete, collect the supernatant that contains gDNA from the elution plate.
6. Store the gDNA sample at 4°C for short-term storage (up to 24 hours), or at -20 or -80°C for long-term storage.